

Name	Model type	Serial number	Date of arrival
Service areas	Yes/No	Test area Check list	OK. Yes/No
Tachometer clean	Y	Tachometer	Y
Imbalance adjustment	Y	Imbalance	Y
Lid latch grease	Y	Lid	Y
Lid springs grease	Y	Rotor recognition	Y
Change Lid seal	Y		
Change Under seal	Y		
Bowl check & clean	Y		
Body clean	Y		
Overlay check & clean	Y		
Rotor areas		Mandatory tests	Pass?
Rotor Check for cracks	Y	PAT Test	Y
Rotor lid seal checks	Y		
Or replace	N		
Lid check for cracks	Y		
Correct number of magnets	Y		
REFRIGERATED MODEL		Mandatory tests	Pass?
Clean Condensor	Y	-10 to 40C	Y
Check Temperature	Y		
Using copy sensor	N		
Check copper lines	Y		
Check wiring	Y		

Further parts or labour required