

Service parts Lists

Year 2017 onward

Ambient & Refrigerated



CENTURION
SCIENTIFIC LIMITED

Why does my Centrifuge need a service?

A yearly service validates your warranty. As with every appliance you purchase servicing extends the lifetime of the product and negates the need for replacement. Which in turn equals money saving

We appreciate how hard your Centrifuge is used and how it is an imperative tool for your Laboratory.

Downtime is therefore not acceptable.

However if it is regularly maintained you improve its reliability and extend its life.

Our Pro-Care parts offer you all that is needed to service our Centrifuges correctly. Full instructions on replacement and any necessary adjustments are included.

Centurion offer this service , to extend your Centrifuge life time and reliability.

As the manufacturer you have , within the package, all the correct manufacturers parts. No waiting and with the correct Instructions.

Using a third party Service company without training and the correct parts is not acceptable. All our distributors for products sold outside of the UK have received training and authority to service our centrifuges.

NOTE.

Centrifuges are extremely dangerous machines.

Although Centurions reputation regarding safety is exceptional.

If a Rotor lid or poor seals are not replaced regularly and they break the damage caused could result in a full unit replacement

Servicing of our centrifuges are also a requirement to validate your warranty and can only be carried out by authorised distributors and using our parts sets only.

Service parts. Ambient, yearly and 5th

Parts sets Pricing

Ambient Models

PrO-	Yearly parts for micro Centrifuge	£49.97
------	-----------------------------------	--------

Care.SP91	5th year parts for micro Centrifuge	£79.97
-----------	-------------------------------------	--------

PrO-

PrO-	Yearly parts for small Centrifuge	£59.97
------	-----------------------------------	--------

Care.SP71	5th year parts for small Centrifuge	£89.97
-----------	-------------------------------------	--------

PrO-

PrO-	Yearly parts for medium Centrifuge	£69.97
------	------------------------------------	--------

Care.SP51	5th year parts for medium Centrifuge	£99.97
-----------	--------------------------------------	--------

PrO-

PrO-Care.SP41	Yearly parts for 2L Centrifuge	£79.97
---------------	--------------------------------	--------

PrO-Care.SP45	5th year parts for 2L Centrifuge	£109.97
---------------	----------------------------------	---------

PrO-Care.SP31	Yearly parts for 3L Centrifuge	£89.97
---------------	--------------------------------	--------

PrO-Care.SP35	5th year parts for 3L Centrifuge	£119.97
---------------	----------------------------------	---------

PrO-Care.SPF1	Yearly parts for Floor Centrifuge	£119.97
---------------	-----------------------------------	---------

PrO-Care.SPF5	5th year parts for Floor Centrifuge	£149.97
---------------	-------------------------------------	---------

Yearly parts as below plus Full Instructions of fitment, adjustment & cleaning

Air spray & Grease spray

Bowl Seal

Lid under Seal

Fuses

5th year parts also include:

Anti vibration mounts, motor boot, overlay and lid springs.

NOTE: If your centrifuge is over 7 years old or old, parts may not be available. We do offer a discount for a replacement centrifuge and we dispose of your old unit as per the Wee directive.

Service parts.Refrigerated, yearly and 5th

Parts sets Pricing

Refrigerated Models

PrO-Care.SPR91	Yearly parts for micro Centrifuge	£69.97
----------------	-----------------------------------	--------

PrO-Care.SPR95	5th year parts for micro Centrifuge	£99.97
----------------	-------------------------------------	--------

PrO-Care.SPR71	Yearly parts for small Centrifuge	£79.97
----------------	-----------------------------------	--------

PrO-Care.SPR75	5th year parts for small Centrifuge	£109.97
----------------	-------------------------------------	---------

PrO-Care.SPR51	Yearly parts for medium Centrifuge	£89.97
----------------	------------------------------------	--------

PrO-Care.SPR55	5th year parts for medium Centrifuge	£129.97
----------------	--------------------------------------	---------

PRO-Care.SPR41	Yearly parts for 2L Centrifuge	£99.97
----------------	--------------------------------	--------

PrO-Care.SPR45	5th year parts for 2L Centrifuge	£129.97
----------------	----------------------------------	---------

PrO-Care.SPR31	Yearly parts for 3L Centrifuge	£109.97
----------------	--------------------------------	---------

PrO-Care.SPR35	5th year parts for 3L Centrifuge	£139.97
----------------	----------------------------------	---------

PrO-Care.SPRF1	Yearly parts for Floor Centrifuge	£139.97
----------------	-----------------------------------	---------

PrO-Care.SPRF5	5th year parts for Floor Centrifuge	£169.97
----------------	-------------------------------------	---------

Yearly parts as below plus Full Instructions of fitment, adjustment & cleaning

Air spray & Grease sparay

Bowl Seal

Under Seal

Fuses

5th year parts also include:

Anti vibration mounts, motor boot, overlay , lid springs & sensor

NOTE: If your centrifuge is over 7 years old or old, parts may not be available. We do offer a discount for a replacement centrifuge and we dispose of your old unit as per the Wee directive.

Service parts. Rotor seals, yearly

Parts sets Pricing

Rotor Seals

PrO-Care	Yearly Replacement Parts	All seals - both inner & outer
RSM1	For Rotors BRK5424, 5436, 5448, 5494	£19.97
RSM2	For Rotors BRK5463 & 5464	£19.97
RSM3	For Rotors BRK5648, 5649, 5698, 5248, 5249 & 5298	£19.97
RSF1	For Rotors BRK5404, 5407, 5408, 5204, 5202	£19.97
RSF2	For Rotors BRK5406, 5412, 5206, 5212	£24.97
RSF3	For Rotors BRK5308, 5324, 5100, 5208, 5224 & 5210	£29.97

PrO-Care	Yearly Replacement Seals	For Sealed Lids
RSSO1	For Rotors BRK1000, (B5250 bucket) & BRK5505	£29.97
RSSO2	For Rotor BRK3000, (B2000, 3000 & 4000 bucket)	£39.97

PrO-Care	1 Year Pack	Important - Replace every 2 months
8011	For Rotor BRK501. Rim gasket pack of 6	£24.98

PrO-Care	Yearly Replacement Seals	All seals - both inner & outer
RSCW	For Rotor 5450. Inner & outer seals. Both tube sizes	£149.97
RSCYT	For Rotors 4420, 4430, 4460. Inner & outer seals	£59.97

NOTE: If your rotor is over 7 years old, it should be replaced.

We do offer a discount for replacement rotors and we offer a disposal service for your old rotor.

PrO-Care

Service, Warranty & Maintenance packages

The lifespan of a centrifuge depends upon how well it is looked after and maintained.

Many of our centrifuges have lived a life of over 15 years.

This is not only due to our stringent manufacturing process and high quality components but also to a yearly servicing package. Making them a cost effective solution of centrifugation.

Whilst our centrifuges come with a confidence giving warranty, yearly services are imperative to ensure your warranty is valid, this also ensures your machine is working to its best ability and lengthens its life span. Also as centrifuges are a dangerous machine, safety checks are required too.

For UK customers we offer a back to base service package including carriage to us and return all within 48 hours. Centrifuge seals , fuses & other parts are replaced, adjustments made, rotor checked for stress cracks, seals changed if necessary (extra charge) PAT (portable appliance test) is included Obligatory in UK.

For worldwide customers we offer a full parts package per model. The designated country distributor services the units accordingly.

Call our Office to discuss or visit our Service website

WWW.ServicePrO-Care.com

www.Centurionscientificglobal.com

Switchboard. 0044 (0) 2392 631225

Research & service

The Old equine Barn
Church Farm
Stoughton Chichester
PO18 9JL. UK

Ambient manufacturing

The Old Grain Barn
Church Farm
Stoughton Chichester
PO18 9JL. UK

Refrigeration Manufacturing

The Old Drier Barn
Church Farm
Stoughton
Chichester
PO18 9JL. UK

Sales

The Old Stables Church
Farm Stoughton
Chichester