

PrO-Care

Yearly service prices. 2020 models onwards



CENTURION
SCIENTIFIC LIMITED

Why does my Centrifuge need a service?

A yearly service validates your warranty

As with every appliance you purchase servicing extends its lifetime and negates the need for constant replacement.

Further we appreciate how hard your centrifuge is used and how it is an imperative tool for your laboratory, downtime is therefore not acceptable. However if it is regularly maintained you improve its reliability and life

Our UK Pro-Care service offers a full service from less than 40 pence a day, a turnaround of 48 hours and as a bonus a free PAT test.

NOTE

Centrifuges can be extremely dangerous machines.

Whilst Centurion's safety record is exceptional, if a seal fails or rotor damage is not checked the damage caused could mean a full replacement.

Using Centurion means the centrifuge utilise's the correct manufacturers parts, our staff have the correct knowledge therefore are fast to service and the charge is lower.

Using a third party without Centurions approval, means incorrect parts fitted, if at all, and invalidates your warranty.

Service prices back to base UK only

Ambient

Micro	£147.00
Small	£247.00
Medium	£267.00
Large2L	£347.00
Large3L	£347.00
Floor	£397.00

Refrigerated

Micro R	£247.00
Small R	£297.00
Medium R	£347.00
Large2L R	£447.00
Large3L R	£447.00
Floor R	£497.00

Prices are inclusive of carriage out and return, a full service including parts, a PAT test, all for as little as 40 pence per day
Prices are for UK mainland only.

Centrifuge and rotor must be correctly packaged individually.
Centurion will not be held responsible for any damage caused by couriers

Turnaround is 48 hours from receipt of unit.

Service includes parts

New seals to relevant areas , fuses, cleaning of tachometer, adjustment of imbalance detector and full individual sensor test.
PAT testing is included for free.

PrO-Care

Service, Warranty & Maintenance packages

The lifespan of a centrifuge depends upon how well it is looked after and maintained.

Many of our centrifuges have lived a life of over 15 years.

This is not only due to our stringent manufacturing process and high quality components but also to a yearly servicing package. Making them a cost effective solution of centrifugation.

Whilst our centrifuges come with a confidence giving warranty, yearly services are imperative to ensure your warranty is valid, this also ensures your machine is working to its best ability and lengthens its life span. Also as centrifuges are a dangerous machine, safety checks are required too.

For UK customers we offer a back to base service package including carriage to us and return all within 48 hours. Centrifuge seals , fuses & other parts are replaced, adjustments made, rotor checked for stress cracks, seals changed if necessary (extra charge) PAT (portable appliance test) is included Obligatory in UK.

For worldwide customers we offer a full parts package per model. The designated country distributor services the units accordingly.

Call our Office to discuss or visit our Service website

WWW.ServicePrO-Care.com

www.Centurionscientificglobal.com

Switchboard. 0044 (0) 2392 631225

Research & service

The Old equine Barn
Church Farm
Stoughton Chichester
PO18 9JL. UK

Ambient manufacturing

The Old Grain Barn
Church Farm
Stoughton Chichester
PO18 9JL. UK

Refrigeration Manufacturing

The Old Drier Barn
Church Farm
Stoughton
Chichester
PO18 9JL. UK

Sales

The Old Stables Church
Farm Stoughton
Chichester